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```
            return;
        }

        $realPath = $value->getRealPath();
        if (!$realPath || !file_exists($realPath) || filesize($realPath) === 0) {
            $fail('The uploaded file is empty or cannot be read.');
```

```
            return;
        }

        // Verify image header and dimensions
        $imageInfo = @getimagesize($realPath);
        if ($imageInfo === false) {
            $fail('The file is not a valid decodable image.');
```

```
            return;
        }

        // Verify supported MIME
        $mime = $imageInfo['mime'] ?? '';
        if (!in_array($mime, ['image/jpeg', 'image/png'], true)) {
            $fail('Only standard JPEG and PNG image formats are supported by the AI engine.');
```

```
            return;
        }

        // Attempt full decode to catch truncated/corrupted payloads
        if (function_exists('imagecreatefromstring')) {
            $content = @file_get_contents($realPath);
            $gd = @imagecreatefromstring($content);
            if ($gd === false) {
                $fail('The image file is corrupted or has an invalid byte stream.');
```

```
                return;
            }
            imagedestroy($gd);
        }
    },
],
'patientId' => ['nullable', 'string', 'max:50'],
'patientName' => ['nullable', 'string', 'max:255'],
'patientAge' => ['nullable', 'integer', 'min:0', 'max:130'],
'patientGender' => ['nullable', 'string', 'max:50'],
'location' => ['nullable', 'string', 'max:100'],
'clinicianNotes' => ['nullable', 'string', 'max:3000'],
];
}

/**
```

```
* Get custom error messages.
*/
public function messages(): array
{
    return [
        'image.required' => 'A kidney scan image is required.',
        'image.file' => 'The uploaded item must be a valid file.',
        'image.image' => 'The file must be a valid image.',
        'image.mimes' => 'Only JPG, JPEG, and PNG images are supported.',
        'image.mimetypes' => 'Only JPG, JPEG, and PNG MIME types (image/jpeg, image/png) are accepted.',
        'image.min' => 'The uploaded file is empty. Please provide a valid ultrasound scan.',
        'image.max' => 'Image size cannot exceed 5MB (model limit).',
        'image.dimensions' => 'Scan dimensions must be between 32x32 and 4096x4096 pixels.',
    ];
}
```

Code Source Intégral 3.2 — Règles de validation (*app/Http/Requests/StoreAnalysisRequest.php*).

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```
*/
public function show(Request $request, int $id): JsonResponse
{
    $analysis = $this->analysisService->getAnalysis(
        $id,
        (int) $request->user()->id,
    );

    $this->auditLogService->log(
        action: 'analysis_viewed',
        resourceType: 'Analysis',
        resourceId: (string) $analysis->id,
        user: $request->user(),
        request: $request
    );

    return response()->json([
        'success' => true,
        'data' => new AnalysisResource($analysis),
    ]);
}

/**
 * Delete an analysis (soft delete).
 *
 * DELETE /api/analyses/{id}
 */
public function destroy(Request $request, int $id): JsonResponse
{
    $this->analysisService->deleteAnalysis(
        $id,
        (int) $request->user()->id,
    );

    $this->auditLogService->log(
        action: 'analysis_deleted',
        resourceType: 'Analysis',
        resourceId: (string) $id,
        user: $request->user(),
        request: $request
    );

    return response()->json([
        'success' => true,
        'message' => 'Analysis deleted successfully.',
    ]);
}

/**
 * Update an analysis (patient metadata / notes).
 *
 * PUT/PATCH /api/analyses/{id}
 */
public function update(\App\Http\Requests\UpdateAnalysisRequest $request, int $id): JsonResponse
{
    $analysis = $this->analysisService->updateAnalysis(
        $id,
        (int) $request->user()->id,
        $request->validated(),
    );
}
```

```

        $this->auditLogService->log(
            action: 'clinician_notes_updated',
            resourceType: 'Analysis',
            resourceId: (string) $analysis->id,
            metadata: [
                'has_notes' => !empty($analysis->clinician_notes),
            ],
            user: $request->user(),
            request: $request
        );

        return response()->json([
            'success' => true,
            'message' => 'Analysis updated successfully.',
            'data' => new AnalysisResource($analysis),
        ]);
    }

    /**
     * Securely stream the analysis ultrasound scan image for authorized owner only.
     *
     * GET /api/analyses/{id}/image
     */
    public function image(Request $request, int $id)
    {
        $analysis = $this->analysisService->getAnalysis(
            $id,
            (int) $request->user()->id,
        );

        $cleanPath = str_replace(['..', '\\'], ['', '/'], (string) $analysis->image_path);
        $path = storage_path('app/public/' . $cleanPath);

        if (!file_exists($path) || is_dir($path)) {
            return response()->json([
                'success' => false,
                'message' => 'Image file not found on server.',
            ], 404);
        }

        return response()->file($path, [
            'Cache-Control' => 'private, max-age=3600',
            'Access-Control-Allow-Origin' => '*',
        ]);
    }

    /**
     * Sign and confirm an analysis by clinician.
     *
     * POST /api/analyses/{id}/confirm
     */
    public function confirm(Request $request, int $id): JsonResponse
    {
        $request->validate([
            'notes' => ['nullable', 'string', 'max:3000'],
        ]);

        $analysis = $this->analysisService->confirmAnalysis(
            $id,
            (int) $request->user()->id,
            $request->input('notes')
        );
    }

```

```

    );

    $this->auditLogService->log(
        action: 'clinician_confirmed_finding',
        resourceType: 'Analysis',
        resourceId: (string) $analysis->id,
        metadata: [
            'confirmed_by' => $request->user()->id,
            'has_notes' => !empty($request->input('notes')),
        ],
        user: $request->user(),
        request: $request
    );

    return response()->json([
        'success' => true,
        'message' => 'Analysis confirmed and signed by clinician.',
        'data' => new AnalysisResource($analysis),
    ]);
}

/**
 * Flag an analysis for specialist triage consultation.
 *
 * POST /api/analyses/{id}/flag
 */
public function flag(Request $request, int $id): JsonResponse
{
    $request->validate([
        'department' => ['nullable', 'string', 'max:50'],
        'reason' => ['nullable', 'string', 'max:1000'],
    ]);

    $department = $request->input('department', 'urology');
    $reason = $request->input('reason');

    $analysis = $this->analysisService->flagAnalysis(
        $id,
        (int) $request->user()->id,
        $department,
        $reason
    );

    $this->auditLogService->log(
        action: 'analysis_flagged_for_consultation',
        resourceType: 'Analysis',
        resourceId: (string) $analysis->id,
        metadata: [
            'department' => $department,
            'reason' => $reason,
        ],
        user: $request->user(),
        request: $request
    );

    return response()->json([
        'success' => true,
        'message' => "Analysis successfully flagged for {$department} triage consultation.",
        'data' => new AnalysisResource($analysis),
    ]);
}

```

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```

        'pathology' => $prediction,
        'class_probabilities' => $analysis->ai_response_payload['class_probabilities'] ?? null,
        'clinical_assessment' => $analysis->ai_response_payload['clinical_assessment'] ?? null,
    ],
    });
}

/**
 * Build a human-readable summary.
 */
private function buildSummary(string $prediction, float $confidence, bool $isUncertain = false): string
{
    $finding = match (strtolower($prediction)) {
        'stone' => 'Nephrolithiasis (Kidney Stone)',
        'cyst' => 'Renal Cyst (Bosniak Classification)',
        'tumor' => 'Suspicious Renal Tissue Mass / Neoplasm',
        default => 'Normal Renal Parenchyma',
    };

    return "AI Multi-Pathology Triage: {$finding} detected with {$confidence}% Detection Confidence. "
        . "This is an automated preliminary screening assessment to assist clinical triage and must be
        confirmed by a qualified medical professional.";
    }

/**
 * Build recommendations based on the prediction.
 */
private function buildRecommendations(string $prediction, float $confidence, bool $isUncertain = false):
array
{
    $base = [
        'All AI triage findings must be reviewed and confirmed by a certified radiologist, nephrologist,
        or urologist.',
        'Correlate findings with clinical history, serum creatinine, GFR, and urinalysis.',
    ];

    return match (strtolower($prediction)) {
        'stone' => array_merge([
            'High suspicion of nephrolithiasis: Urological consultation recommended.',
            'Assess stone dimensions, location, and potential urinary tract obstruction.',
            'Consider non-contrast CT abdomen/pelvis if ultrasound was initial modality.',
            'Ensure adequate patient hydration and pain management per clinical protocol.',
        ], $base),
        'cyst' => array_merge([
            'Cortical renal cyst identified consistent with benign fluid-filled architecture (Bosniak
            I/II).',
            'Periodic annual ultrasound monitoring recommended to verify dimensional stability.',
            'Correlate with renal function panel; surgical intervention typically unnecessary unless
            symptomatic.',
        ], $base),
        'tumor' => array_merge([
            'CRITICAL ALERT: Solid heterogeneous renal mass detected. Immediate specialist referral
            required.',
            'Order urgent contrast-enhanced multiphasic CT abdomen/pelvis or renal MRI with contrast.',
            'Schedule urgent urological oncology consultation for staging and histological
            characterization.',
        ], $base),
        'normal' => array_merge([
            'Renal parenchyma appearance is within normal limits; no nephrolithiasis, cyst, or solid
            mass detected.',
            'If clinical symptoms persist, consider differential diagnosis and non-contrast CT.',
        ], $base);
    };
}

```

```
        'Maintain routine preventative renal care and hydration.',  
        ], $base),  
        default => $base,  
    };  
}  
}
```

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