DOCUMENT TECHNIQUE N° 4 — FRONTEND CLIENT REACT 18 & TYPESCRIPT

KIDNEYVISION AI — INTERFACE CLINIQUE MODERNE

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2. Typage Statique Strict & Interfaces TypeScript

```
/**
 * Central TypeScript Type Definitions Hub
 */

export * from './auth';
export * from './analysis';
export * from './dashboard';
export * from './settings';
export * from './health';
export * from './scan';
export * from './analytics';
```

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```
import { useState } from "react";
import { BrowserRouter, Routes, Route, Navigate, useLocation, Link, useNavigate } from "react-router-dom";
import {
  Activity, ShieldCheck, History, Settings, LogOut, FileText,
  LayoutDashboard, User, HelpCircle, Eye, Menu, X, PlusCircle, BrainCircuit
} from "lucide-react";

import LandingPage from "../components/LandingPage";
import LoginPage from "../components/auth/LoginPage";
import RegisterPage from "../components/auth/RegisterPage";
import ProtectedRoute from "../components/auth/ProtectedRoute";

import React, { Suspense } from "react";
import LoadingSkeleton from "../components/common/LoadingSkeleton";
import ErrorBoundary from "../components/common/ErrorBoundary";

const AppDashboard = React.lazy(() => import("../components/AppDashboard"));
const AppNewAnalysis = React.lazy(() => import("../components/AppNewAnalysis"));
const AppHistory = React.lazy(() => import("../components/AppHistory"));
const AppReports = React.lazy(() => import("../components/AppReports"));
const AppSettings = React.lazy(() => import("../components/AppSettings"));
const AppHelpCenter = React.lazy(() => import("../components/AppHelpCenter"));
const GuestAnalysis = React.lazy(() => import("../components/GuestAnalysis"));
const ForgotPassword = React.lazy(() => import("../components/ForgotPassword"));
const ResetPassword = React.lazy(() => import("../components/ResetPassword"));
const PrivacyPage = React.lazy(() => import("../components/PrivacyPage"));
const TermsPage = React.lazy(() => import("../components/TermsPage"));

import { AuthProvider, useAuth } from "../contexts/AuthContext";

function PortalLayout() {
  const { user, logout } = useAuth();
  const [mobileMenuOpen, setMobileMenuOpen] = useState(false);
  const location = useLocation();
  const navigate = useNavigate();

  const handleLogout = () => {
    logout();
  };

  const getActiveTab = () => {
    const path = location.pathname;
    if (path.includes("/analysis")) return "analysis";
    if (path.includes("/history")) return "history";
    if (path.includes("/reports")) return "reports";
    if (path.includes("/settings")) return "settings";
    if (path.includes("/help")) return "help";
    return "dashboard";
  };

  const activeTab = getActiveTab();

  return (
```

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
5. Module de Diagnostic & Glisser-Déposer (AppNewAnalysis.tsx)

```

/**
 * @license
 * SPDX-License-Identifier: Apache-2.0
 */

import React, { useState, useRef } from "react";
import {
  Upload, Sparkles, Loader2, Play, CircleDot, ShieldCheck,
  Eye, FileText, CheckCircle2, ChevronRight, Download, Dna,
  ArrowRight, Heart, AlertCircle, FileDown, X, RefreshCw,
  User, Calendar, Layers, Image as ImageIcon, Check, Sliders,
  HelpCircle, Info, Stethoscope, Hash, Trash2, Crosshair,
  Ruler, Activity, Cpu, CheckCircle, Maximize2
} from "lucide-react";
import { createPrediction, downloadReport, updateAnalysis, getReportPdfBlobUrl } from "../services";
import { DiagnosisResult, KidneyAnalysis } from "../types";
import { validateScanFile } from "../utils/fileValidation";
import { parseClinicalError, ClinicalError } from "../utils/errorParser";

interface AppNewAnalysisProps {
  onAddSuccess: () => void;
}

export default function AppNewAnalysis({ onAddSuccess }: AppNewAnalysisProps) {
  // Input fields state – starts completely empty for clinical doctor entry
  const [patientId, setPatientId] = useState("");
  const [patientName, setPatientName] = useState("");
  const [patientAge, setPatientAge] = useState("");
  const [patientGender, setPatientGender] = useState<"Male" | "Female" | "Other">("Male");
  const [location, setLocation] = useState<"Left Kidney" | "Right Kidney" | "Unspecified">("Left Kidney");
  const [clinicianNotes, setClinicianNotes] = useState("");

  // Custom file upload state
  const [uploadedFile, setUploadedFile] = useState(null);
  const [uploadedFilePreview, setUploadedFilePreview] = useState(null);
  const [uploadError, setUploadError] = useState(null);
  const [isDragging, setIsDragging] = useState(false);
  const fileInputRef = useRef(null);

  // View mode for results: "panel" (Dual Clinical Panel) | "gradcam" | "cleaned" | "raw"
  const [heatmapViewMode, setHeatmapViewMode] = useState<"panel" | "gradcam" | "cleaned" | "raw">("panel");

  // Diagnosis & Pipeline states
  const [isScanning, setIsScanning] = useState(false);
  const [currentStageIdx, setCurrentStageIdx] = useState(0);
  const [clinicalError, setClinicalError] = useState(null);
  const [activeAnalysis, setActiveAnalysis] = useState(null);

  // Clinician notes save state
  const [isSavingNotes, setIsSavingNotes] = useState(false);
  const [noteSavedSuccess, setNoteSavedSuccess] = useState(false);
  const [pdfPreviewUrl, setPdfPreviewUrl] = useState(null);
  const [isPreviewLoading, setIsPreviewLoading] = useState(false);

  // Analysis Pipeline Stages
  const ANALYSIS_STAGES = [
    { id: "uploading", label: "Uploading Scan", desc: "Streaming acoustic tensor to secure medical

```

```

pipeline..." },
  { id: "gatekeeper", label: "Stage-1 Gatekeeper Screening", desc: "MobileNetV3 evaluating renal vs non-
kidney scan..." },
  { id: "preprocessing", label: "Spatial Cleaning & Inpainting", desc: "Top-Hat filtering & Telea
inpainting on 224x224 tensor..." },
  { id: "analyzing", label: "ConvNext 4-Class Inference", desc: "Evaluating Cyst, Normal, Stone, Tumor
pathologies..." },
  { id: "generating", label: "Grad-CAM Heatmap & Biometrics", desc: "Synthesizing spatial feature
activation map & lesion diameter..." },
  { id: "completed", label: "Diagnostic Report Ready", desc: "Audit logged and ready for radiologist
confirmation." },
];

// Process chosen file
const processSelectedFile = async (file: File) => {
  setUploadError(null);
  const validation = await validateScanFile(file);
  if (!validation.valid) {
    setUploadError(validation.error || "The selected file is not a valid medical scan image.");
    setUploadedFile(null);
    setUploadedFilePreview(null);
    if (fileInputRef.current) fileInputRef.current.value = "";
    return;
  }

  setUploadedFile(file);
  const url = URL.createObjectURL(file);
  setUploadedFilePreview(url);
  setActiveAnalysis(null);

  // If patient name is currently blank, suggest an automated placeholder from filename
  if (!patientName.trim()) {
    const cleanName = file.name.replace(/\.([^/]+)$/g, "").replace(/[-_]/g, " ");
    // Only set if not a generic generic hash
    if (cleanName.length > 2 && !cleanName.match(/^[0-9a-f]{8,}$/i)) {
      setPatientName(cleanName.charAt(0).toUpperCase() + cleanName.slice(1));
    }
  }
};

const handleImageFileChange = async (e: React.ChangeEvent) => {
  if (e.target.files && e.target.files[0]) {
    await processSelectedFile(e.target.files[0]);
  }
};

const handleDragOver = (e: React.DragEvent) => {
  e.preventDefault();
  e.stopPropagation();
  setIsDragging(true);
};

const handleDragLeave = (e: React.DragEvent) => {
  e.preventDefault();
  e.stopPropagation();
  setIsDragging(false);
};

const handleDrop = async (e: React.DragEvent) => {
  e.preventDefault();
  e.stopPropagation();

```

```

    setIsDragging(false);
    if (e.dataTransfer.files && e.dataTransfer.files[0]) {
      await processSelectedFile(e.dataTransfer.files[0]);
    }
  };

const handleRemoveFile = () => {
  setUploadedFile(null);
  if (uploadedFilePreview) {
    URL.revokeObjectURL(uploadedFilePreview);
    setUploadedFilePreview(null);
  }
  setUploadError(null);
  setActiveAnalysis(null);
  if (fileInputRef.current) fileInputRef.current.value = "";
};

// Run AI Scan Execution
const handleTriggerAnalysis = async () => {
  if (!uploadedFile && !uploadedFilePreview) {
    setUploadError("Please upload a kidney ultrasound scan before running analysis.");
    return;
  }

  setIsScanning(true);
  setCurrentStageIdx(0);
  setClinicalError(null);
  setActiveAnalysis(null);

  const stageTimer = setInterval(() => {
    setCurrentStageIdx((prev) => {
      if (prev < 3) return prev + 1;
      return prev;
    });
  }, 600);

  try {
    const finalPatientName = patientName.trim() || "Unspecified Patient";
    const finalPatientId = patientId.trim() || undefined;
    const parsedAge = patientAge.trim() ? parseInt(patientAge, 10) : undefined;

    const res = await createPrediction(uploadedFile, {
      patientId: finalPatientId,
      patientName: finalPatientName,
      patientAge: parsedAge,
      patientGender,
      location,
      imageUrl: uploadedFilePreview || undefined,
      clinicianNotes: clinicianNotes.trim() || undefined,
    });

    clearInterval(stageTimer);
    setCurrentStageIdx(4);

    setTimeout(() => {
      setActiveAnalysis(res);
      if (res.clinicianNotes) {
        setClinicianNotes(res.clinicianNotes);
      }
      setIsScanning(false);
      const isNormal = (res.prediction || '').toLowerCase() === 'normal' || res.diagnosis ===

```

```

DiagnosisResult.NORMAL_FINDINGS;
    setHeatmapViewMode(isNormal ? (res.cleanedScanUrl ? "cleaned" : "raw") : (res.visualizationUrl ?
"panel" : res.heatmapUrl ? "gradcam" : "raw"));
    }, 400);
  } catch (e: any) {
    clearInterval(stageTimer);
    setIsScanning(false);
    const parsed = e?.clinicalError || parseClinicalError(e);
    setClinicalError(parsed);
  }
};

// Save report action (persists clinician notes and metadata updates)
const handleSaveReport = async () => {
  if (!activeAnalysis?.id) return;
  setIsSavingNotes(true);
  try {
    await updateAnalysis(activeAnalysis.id, {
      clinicianNotes: clinicianNotes.trim()
    });
    setNoteSavedSuccess(true);
    setTimeout(() => setNoteSavedSuccess(false), 3000);
  } catch (e: any) {
    alert("Error saving clinician notes: " + e.message);
  } finally {
    setIsSavingNotes(false);
  }
};

return (

  /* Top Header & Breadcrumbs matching Dashboard style */

```

AI Computer Vision Inference Console

New Renal Diagnostic Analysis

Upload renal CT or ultrasound scans to run Stage-1 Gatekeeper screening, ConvNeXt-Tiny 4-class classification (Cyst, Normal, Stone, Tumor), and Grad-CAM explainability localization.

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```
import React, { useState, useMemo, useEffect } from 'react';
import { useNavigate } from 'react-router-dom';
import { DashboardSidebar } from '../components/dashboard/DashboardSidebar';
import { DashboardHeader } from '../components/dashboard/DashboardHeader';
import { KpiGrid } from '../components/dashboard/KpiGrid';
import { Kidney3DViewer } from '../components/dashboard/Kidney3DViewer';
import { DiagnosticOverview } from '../components/dashboard/DiagnosticOverview';
import { ScanVolumeChart } from '../components/dashboard/ScanVolumeChart';
import { ConfidenceChart } from '../components/dashboard/ConfidenceChart';
import { RiskDistribution } from '../components/dashboard/RiskDistribution';
import { LateralityChart } from '../components/dashboard/LateralityChart';
import { DemographicsChart } from '../components/dashboard/DemographicsChart';
import { RecentClinicalActivity } from '../components/dashboard/RecentClinicalActivity';
import { ScanDetailModal } from '../components/dashboard/ScanDetailModal';

import {
  computeDashboardMetrics,
  computeDiagnosticDistribution,
  computeRiskDistribution,
  computeLateralityDistribution,
  computeConfidenceBuckets,
  computeTimelineData,
  computeAgeCohorts,
  filterScansByDateRange,
} from '../data/mockAnalytics';

import { getAnalyses } from '../services/analysis.service';
import { KidneyAnalysis, DiagnosisResult } from '../types/analysis';

import type { DateRangeOption } from '../types/dashboard';
import type { Diagnosis, ScanRecord, ScanStatus } from '../types/scan';
import { MOCK_SCANS } from '../data/mockScans';
import { Menu, Upload, CheckCircle2, ArrowRight } from 'lucide-react';

export interface DashboardProps {
  hideSidebar?: boolean;
}

function mapAnalysisToScanRecord(analysis: KidneyAnalysis): ScanRecord {
  const confidence = typeof analysis.confidence === 'number' ? Number(analysis.confidence.toFixed(1)) :
  null;
  const findingLower = (analysis.pathologyFinding || '').toLowerCase();
  const predLower = (analysis.prediction || '').toLowerCase();
  const isRejected = analysis.diagnosis === DiagnosisResult.REJECTED || findingLower.includes('rejected') ||
  findingLower.includes('invalid') || findingLower.includes('non-kidney');
  const isNormal = predLower === 'normal' || findingLower.includes('normal') || analysis.diagnosis ===
  DiagnosisResult.NORMAL_FINDINGS;
  const isPending = !isNormal && (analysis.diagnosis === DiagnosisResult.ANALYZING || analysis.diagnosis ===
  DiagnosisResult.REVIEW_REQUIRED);

  let diagnosis: Diagnosis = 'normal';
  let riskLevel: 'low' | 'elevated' | 'urgent' = 'low';

  if (isRejected) {
    diagnosis = 'rejected';
    riskLevel = 'elevated';
  } else if (isNormal) {
```

```

    diagnosis = 'normal';
    riskLevel = 'low';
  } else if (findingLower.includes('tumor') || predLower === 'tumor') {
    diagnosis = 'tumor';
    riskLevel = 'urgent';
  } else if (findingLower.includes('stone') || predLower === 'stone' || analysis.diagnosis ===
DiagnosisResult.ANOMALY_DETECTED) {
    diagnosis = 'nephrolithiasis';
    riskLevel = confidence && confidence > 90 ? 'urgent' : 'elevated';
  } else if (findingLower.includes('cyst') || predLower === 'cyst') {
    diagnosis = 'cyst';
    riskLevel = 'low';
  } else if (isPending) {
    diagnosis = 'pending';
    riskLevel = 'elevated';
  }
}

const rawId = String(analysis.id || '');
const id = rawId.startsWith('KV-') || rawId.startsWith('PT-') ? rawId : `KV-${rawId.padStart(3, '0')}`;

let formattedDate = 'Recent';
let date = new Date().toISOString().slice(0, 10);
if (analysis.createdAt) {
  try {
    const d = new Date(analysis.createdAt);
    if (!isNaN(d.getTime())) {
      date = d.toISOString().slice(0, 10);
      formattedDate = d.toLocaleDateString('en-US', { month: 'short', day: 'numeric', year: 'numeric' });
    }
  } catch {
    // keep defaults
  }
}

const loc = String(analysis.location || '').toLowerCase();
const laterality = loc.includes('right') ? 'right' : loc.includes('bilateral') ? 'bilateral' : 'left';
const diamMm = analysis.clinicalAssessment?.estimated_diameter_mm || analysis.dimensions?.length;

const rawStatus = (analysis.status || '').toLowerCase();
let scanStatus: ScanStatus = 'reviewed';
if (
  rawStatus === 'review' ||
  rawStatus === 'pending' ||
  analysis.diagnosis === DiagnosisResult.REVIEW_REQUIRED ||
  analysis.isDraft
) {
  scanStatus = 'under_review';
} else if (rawStatus === 'flagged') {
  scanStatus = 'flagged';
} else if (rawStatus === 'confirmed' || rawStatus === 'completed') {
  scanStatus = 'reviewed';
} else {
  scanStatus = analysis.confirmedAt ? 'reviewed' : 'under_review';
}

return {
  id,
  patientAlias: analysis.patientName || analysis.patientId || `Patient ${rawId}`,
  date,
  formattedDate,
  laterality,

```

```

    diagnosis,
    confidence,
    riskLevel,
    status: scanStatus,
    patientAge: analysis.patientAge || 45,
    patientSex: String(analysis.patientGender || '').toLowerCase() === 'female' ? 'female' : 'male',
    kidneyPole: 'lower',
    stoneSizeMm: diamMm,
    hasGradCam: !isNormal && Boolean(analysis.heatmapUrl || analysis.visualizationUrl),
    notes: analysis.clinicianNotes || analysis.recommendation || 'Clinical record from database.',
    confirmedAt: analysis.confirmedAt,
    confirmedByName: analysis.confirmedByName,
    flaggedFor: analysis.flaggedFor,
    flaggedAt: analysis.flaggedAt,
    pathologyTitle: analysis.pathologyFinding,
    classProbabilities: analysis.classProbabilities,
    clinicalAssessment: analysis.clinicalAssessment,
    gateEvaluation: analysis.gateEvaluation,
    visualizationUrl: isNormal ? undefined : analysis.visualizationUrl,
    cleanedScanUrl: analysis.cleanedScanUrl,
    imageUrl: analysis.imageUrl,
    heatmapUrl: isNormal ? undefined : analysis.heatmapUrl,
    primaryDiagnosis: analysis.primaryDiagnosis,
    modelVersion: analysis.modelVersion,
    telemetry: analysis.telemetry,
  };
}

export const Dashboard: React.FC = ({ hideSidebar = false }) => {
  const navigate = useNavigate();
  const [isOpenMobile, setIsOpenMobile] = useState(false);

  // Backend analyses state
  const [backendScans, setBackendScans] = useState([]);

  // Filter state
  const [dateRange, setDateRange] = useState('30D');
  const [selectedDiagnosis, setSelectedDiagnosis] = useState('all');

  // Modal inspection state
  const [selectedScan, setSelectedScan] = useState(null);

  // New Analysis Demo State
  const [showNewScanModal, setShowNewScanModal] = useState(false);
  const [newScanStep, setNewScanStep] = useState(1);

  // Load real analyses from Laravel API if available
  const loadBackendScans = async () => {
    try {
      const res = await getAnalyses(1, 100);
      if (res && res.data && res.data.length > 0) {
        const mapped = res.data.map(mapAnalysisToScanRecord);
        setBackendScans(mapped);
      }
    } catch {
      // Backend offline or user not logged in; gracefully proceed with MOCK_SCANS
    }
  };

  useEffect(() => {
    loadBackendScans();
  }, []);

```

```
}, []]);

// Real scans from backend database with graceful fallback to MOCK_SCANS
const allScans = useMemo(() => {
  return backendScans.length > 0 ? backendScans : MOCK_SCANS;
}, [backendScans]);

// Filter scans by active date range
const filteredRangeScans = useMemo(() => {
  return filterScansByDateRange(allScans, dateRange);
}, [allScans, dateRange]);

// Derived metrics & chart datasets
const metrics = useMemo(() => {
  return computeDashboardMetrics(filteredRangeScans);
}, [filteredRangeScans]);

const diagnosticDistribution = useMemo(() => {
  return computeDiagnosticDistribution(filteredRangeScans);
}, [filteredRangeScans]);

const riskDistribution = useMemo(() => {
  return computeRiskDistribution(filteredRangeScans);
}, [filteredRangeScans]);

const lateralityDistribution = useMemo(() => {
  return computeLateralityDistribution(filteredRangeScans);
}, [filteredRangeScans]);

const confidenceBuckets = useMemo(() => {
  return computeConfidenceBuckets(filteredRangeScans);
}, [filteredRangeScans]);

const timelineData = useMemo(() => {
  return computeTimelineData(filteredRangeScans);
}, [filteredRangeScans]);

const ageCohorts = useMemo(() => {
  return computeAgeCohorts(filteredRangeScans);
}, [filteredRangeScans]);

const handleScanUpdated = (updatedScan: ScanRecord) => {
  setBackendScans((prev) =>
    prev.map((s) => (s.id === updatedScan.id ? updatedScan : s))
  );
  setSelectedScan(updatedScan);
};

return (
  <>
    </>
  </>

```

```
    { /* Mobile Top Header Bar (only when standalone) */ }
    { !hideSidebar && (
```

```
                                setIsOpenMobile(true)}
    className="p-2 text-slate-600 hover:text-slate-900 rounded-lg border border-slate-200"
    aria-label="Open sidebar menu"
    >
```

KidneyVision AI

MS

```
    )}
```

```
    { /* 2. Top Header with Controls */ }
    {
        loadBackendScans();
        setDateRange((prev) => prev);
    }
    pendingCount={metrics.pendingReviewCount}
    onNewScanClick={() => setShowNewScanModal(true)}
    />
```

```
    { /* Dashboard Body Container */ }
```


```
    { /* 3. Top KPI Cards */ }
```

```
    { /* 4. Main Visual Section: 3D Kidney & Diagnostic Overview */ }
```

```
    { /* Left 3D Kidney Viewer (compact 4 cols on lg) */ }
```

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```

/**
 * @license
 * SPDX-License-Identifier: Apache-2.0
 */

import React, { useState, useEffect } from "react";
import { Link } from "react-router-dom";
import {
  Search, Calendar, Filter, Eye, Trash2, Sliders, X,
  ChevronRight, ChevronLeft, Loader2, Award, ClipboardCheck, ArrowUpRight, CheckCircle2, AlertCircle,
  RefreshCw,
  FileText, Download, FileDown, Sparkles, ZoomIn, ExternalLink, Crosshair, Ruler, Activity, Cpu, Layers,
  ShieldCheck
} from "lucide-react";
import { getAnalyses, deleteAnalysis, downloadReport, updateAnalysis, getReportPdfBlobUrl } from
"../services";
import { KidneyAnalysis, DiagnosisResult } from "../types";
import { parseClinicalError, ClinicalError } from "../utils/errorParser";
import LoadingSkeleton from "../common/LoadingSkeleton";

export default function AppHistory() {
  const [analyses, setAnalyses] = useState([]);
  const [loading, setLoading] = useState(true);
  const [searchId, setSearchId] = useState("");
  const [statusFilter, setStatusFilter] = useState("All Results");
  const [showFilterDropdown, setShowFilterDropdown] = useState(false);
  const [dateFilter, setDateFilter] = useState<"All Time" | "Last 7 Days" | "Last 30 Days" | "Last 90 Days">
("All Time");
  const [showDateDropdown, setShowDateDropdown] = useState(false);
  const [currentPage, setCurrentPage] = useState(1);
  const [meta, setMeta] = useState({ current_page: 1, last_page: 1, total: 0 });

  const [selectedAnalysis, setSelectedAnalysis] = useState(null);
  const [dossierViewMode, setDossierViewMode] = useState<"panel" | "gradcam" | "cleaned" | "raw">("panel");
  const [submittingNote, setSubmittingNote] = useState(false);
  const [editingNoteText, setEditingNoteText] = useState("");
  const [isEditingNote, setIsEditingNote] = useState(false);
  const [noteSavedFeedback, setNoteSavedFeedback] = useState(false);

  // Notification Banner state
  const [notification, setNotification] = useState<{ type: 'error' | 'success' | 'info'; title: string;
message: string; suggestion?: string } | null>(null);

  // PDF Preview and Download states
  const [previewPdfUrl, setPreviewPdfUrl] = useState(null);
  const [previewPdfTitle, setPreviewPdfTitle] = useState("");
  const [previewLoadingId, setPreviewLoadingId] = useState(null);
  const [downloadLoadingId, setDownloadLoadingId] = useState(null);

  // Enlarged Image Modal state
  const [enlargedImageUrl, setEnlargedImageUrl] = useState(null);
  const [enlargedImageTitle, setEnlargedImageTitle] = useState("");

  const handlePreviewPdf = async (item: KidneyAnalysis, e?: React.MouseEvent) => {
    if (e) e.stopPropagation();
    setPreviewLoadingId(item.id);
    try {

```

```

    const url = await getReportPdfBlobUrl(item.id);
    setPreviewPdfTitle(`Diagnostic Report - ${item.patientName} || item.patientId} (${item.id})`);
    setPreviewPdfUrl(url);
  } catch (err: any) {
    const parsed = parseClinicalError(err);
    setNotification({
      type: 'error',
      title: parsed.title,
      message: parsed.message,
      suggestion: parsed.suggestion
    });
  } finally {
    setPreviewLoadingId(null);
  }
};

```

```

const handleDownloadPdf = async (item: KidneyAnalysis, e?: React.MouseEvent) => {
  if (e) e.stopPropagation();
  setDownloadLoadingId(item.id);
  try {
    await downloadReport(item.id, `kidneyvision_report_${item.id}.pdf`);
    setNotification({
      type: 'success',
      title: 'Report Downloaded',
      message: `Diagnostic PDF report for study #${item.id} downloaded successfully.`
    });
  } catch (err: any) {
    const parsed = parseClinicalError(err);
    setNotification({
      type: 'error',
      title: parsed.title,
      message: parsed.message,
      suggestion: parsed.suggestion
    });
  } finally {
    setDownloadLoadingId(null);
  }
};

```

```

const openDetails = (item: KidneyAnalysis) => {
  setSelectedAnalysis(item);
  setDossierViewMode(item.visualizationUrl ? "panel" : item.heatmapUrl ? "gradcam" : "raw");
  setEditingNoteText(item.clinicianNotes || "");
  setIsEditingNote(false);
  setNoteSavedFeedback(false);
};

```

```

async function loadAnalyses(page: number) {
  setLoading(true);
  try {
    const response = await getAnalyses(page, 10);
    setAnalyses(response.data);
    setMeta(response.meta);
    setCurrentPage(response.meta.current_page);
  } catch (e: any) {
    const parsed = parseClinicalError(e);
    setNotification({
      type: 'error',
      title: parsed.title,
      message: parsed.message,
      suggestion: parsed.suggestion
    });
  }
}

```

```

    });
  } finally {
    setLoading(false);
  }
}

useEffect(() => {
  loadAnalyses(currentPage);
}, [currentPage]);

// Handle Deletion (DELETE /analyses/:id)
const handleDelete = async (id: string, e: React.MouseEvent) => {
  e.stopPropagation(); // Avoid triggering details modal drawer
  if (confirm(`Are you sure you want to delete patient record ${id}?`)) {
    try {
      await deleteAnalysis(id);
      setAnalyses(prev => prev.filter(x => x.id !== id));
      setMeta(prev => ({ ...prev, total: Math.max(0, prev.total - 1) }));
      if (selectedAnalysis?.id === id) {
        setSelectedAnalysis(null);
      }
      setNotification({
        type: 'success',
        title: 'Record Deleted',
        message: `Patient analysis record #${id} was deleted successfully.`
      });
    } catch (err: any) {
      const parsed = parseClinicalError(err);
      setNotification({
        type: 'error',
        title: parsed.title,
        message: parsed.message,
        suggestion: parsed.suggestion
      });
    }
  }
};

// Perform search & state filtering
const filteredAnalyses = analyses.filter((item) => {
  const query = searchId.trim().toLowerCase();
  const matchesSearch = !query ||
    (item.patientId && item.patientId.toLowerCase().includes(query)) ||
    (item.patientName && item.patientName.toLowerCase().includes(query)) ||
    (item.id && String(item.id).toLowerCase().includes(query)) ||
    (item.location && item.location.toLowerCase().includes(query));

  // Status filter - supporting all 4 pathologies and triage states
  let matchesStatus = true;
  const itemDiag = item.diagnosis;
  const itemFinding = (item.pathologyFinding || item.prediction || '').toLowerCase();

  if (statusFilter === "Normal Findings") {
    matchesStatus = itemDiag === DiagnosisResult.NORMAL_FINDINGS || itemFinding.includes('normal');
  } else if (statusFilter === "Kidney Stone") {
    matchesStatus = itemDiag === DiagnosisResult.STONE || itemFinding.includes('stone') ||
    itemFinding.includes('nephro');
  } else if (statusFilter === "Renal Cyst") {
    matchesStatus = itemDiag === DiagnosisResult.CYST || itemFinding.includes('cyst');
  } else if (statusFilter === "Renal Tumor") {
    matchesStatus = itemDiag === DiagnosisResult.TUMOR || itemFinding.includes('tumor');
  }
});

```

```

    } else if (statusFilter === "Review Required") {
        matchesStatus = itemDiag === DiagnosisResult.REVIEW_REQUIRED;
    } else if (statusFilter === "Scan Rejected") {
        matchesStatus = itemDiag === DiagnosisResult.REJECTED || itemFinding.includes('reject') ||
itemFinding.includes('invalid');
    }

// Date filter
let matchesDate = true;
if (dateFilter !== "All Time") {
    const itemTime = new Date(item.createdAt).getTime();
    const now = Date.now();
    if (!isNaN(itemTime)) {
        if (dateFilter === "Last 7 Days") matchesDate = itemTime >= now - 7 * 86400000;
        else if (dateFilter === "Last 30 Days") matchesDate = itemTime >= now - 30 * 86400000;
        else if (dateFilter === "Last 90 Days") matchesDate = itemTime >= now - 90 * 86400000;
    }
}

return matchesSearch && matchesStatus && matchesDate;
});

// Relative date formatting
const formatDate = (isoString: string) => {
    try {
        const date = new Date(isoString);
        if (isNaN(date.getTime())) {
            return { date: "Unknown Date", time: "--:--" };
        }
        const months = ["Jan", "Feb", "Mar", "Apr", "May", "Jun", "Jul", "Aug", "Sep", "Oct", "Nov", "Dec"];
        const formattedDate = `${months[date.getMonth()]} ${date.getDate()}, ${date.getFullYear()}`;

        let hours = date.getHours();
        const minutes = date.getMinutes().toString().padStart(2, '0');
        const ampm = hours >= 12 ? 'PM' : 'AM';
        hours = hours % 12;
        hours = hours ? hours : 12;

        return {
            date: formattedDate,
            time: `${hours}:${minutes} ${ampm}`
        };
    } catch {
        return { date: "Unknown Date", time: "--:--" };
    }
};

return (
    { /* Upper header section */ }

```

Analysis History

Review past ultrasound analyses, diagnostic confidence scores, and historical patient imaging records.

```
{/* Top bar controls: Search & Filters (identical to Screen 3 Header controls) */}
```

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