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```

"""
KidneyVision AI – Stage 1 Gate Model Training Script
Binary AI Filter: Kidney Scan vs. Non-Kidney (Chest X-ray, Bone X-ray, Dental X-ray, etc.)

Usage:
    python train_gate_model.py --dataset-dir /path/to/dataset_gate --output-model kidney_gate_model.keras
"""

import os
import sys
import glob
import hashlib
import random
import argparse
import numpy as np
from PIL import Image

import tensorflow as tf
from tensorflow.keras import layers, models, callbacks, optimizers, metrics

# Reproducibility
SEED = 42
os.environ["PYTHONHASHSEED"] = str(SEED)
random.seed(SEED)
np.random.seed(SEED)
tf.random.set_seed(SEED)

def clean_and_audit_dataset(directory_path: str):
    """Scans directory to remove unreadable, zero-byte, and duplicate images."""
    if not os.path.exists(directory_path):
        print(f"[WARNING] Directory not found: {directory_path}. Skipping audit.")
        return

    valid_extensions = (".jpg", ".jpeg", ".png", ".bmp", ".tif", ".tiff", ".jfif")
    seen_hashes = set()
    removed_corrupted = 0
    removed_duplicates = 0
    total_valid = 0

    print(f"\n[DATA AUDIT] Auditing images in: {directory_path}...")

    for root, _, files in os.walk(directory_path):
        for fname in files:
            fpath = os.path.join(root, fname)

            if not fname.lower().endswith(valid_extensions):
                continue

            # Check file size > 0
            if os.path.getsize(fpath) == 0:
                try:
                    os.remove(fpath)
                    removed_corrupted += 1
                except OSError:
                    pass
                continue

            # Verify image format integrity
            try:
                with Image.open(fpath) as img:

```

```

        img.verify()

        # MD5 Hash check for duplicate removal
        with open(fpath, "rb") as f:
            file_hash = hashlib.md5(f.read()).hexdigest()

        if file_hash in seen_hashes:
            try:
                os.remove(fpath)
                removed_duplicates += 1
            except OSError:
                pass
        else:
            seen_hashes.add(file_hash)
            total_valid += 1

    except Exception:
        try:
            os.remove(fpath)
            removed_corrupted += 1
        except OSError:
            pass

    print(f"Audit Complete for {directory_path}:")
    print(f" - Valid unique images: {total_valid}")
    print(f" - Corrupted files removed: {removed_corrupted}")
    print(f" - Duplicate files removed: {removed_duplicates}")

def build_gate_model(input_shape=(224, 224, 3)):
    """Builds MobileNetV3-Small binary classification model."""
    base_model = tf.keras.applications.MobileNetV3Small(
        input_shape=input_shape,
        include_top=False,
        weights="imagenet"
    )
    base_model.trainable = False # Frozen for Phase 1

    inputs = tf.keras.Input(shape=input_shape, name="scan_input")
    x = tf.keras.applications.mobilenet_v3.preprocess_input(inputs)
    x = base_model(x, training=False)

    # Clinical Triage Head
    x = layers.GlobalAveragePooling2D(name="gap")(x)
    x = layers.BatchNormalization(name="gate_bn1")(x)
    x = layers.Dropout(0.35, name="gate_dropout1")(x)
    x = layers.Dense(128, activation="relu", name="gate_dense1")(x)
    x = layers.BatchNormalization(name="gate_bn2")(x)
    x = layers.Dropout(0.2, name="gate_dropout2")(x)
    outputs = layers.Dense(1, activation="sigmoid", name="kidney_probability")(x)

    model = tf.keras.Model(inputs=inputs, outputs=outputs, name="KidneyVision_Gate_MobileNetV3")
    return model, base_model

def main():
    parser = argparse.ArgumentParser(description="Train KidneyVision AI Gate Model with Hard-Negative X-
    Rays")
    parser.add_argument("--dataset-dir", type=str, default="dataset_gate", help="Path to gate dataset
    directory")
    parser.add_argument("--train-dir", type=str, default=None, help="Explicit path to train directory")

```

```

parser.add_argument("--val-dir", type=str, default=None, help="Explicit path to val directory")
parser.add_argument("--img-size", type=int, default=224, help="Target image dimension (default: 224)")
parser.add_argument("--batch-size", type=int, default=32, help="Batch size (default: 32)")
parser.add_argument("--epochs-p1", type=int, default=8, help="Phase 1 epochs (frozen backbone)")
parser.add_argument("--epochs-p2", type=int, default=12, help="Phase 2 epochs (fine-tuning)")
parser.add_argument("--threshold", type=float, default=0.75, help="Strict admission threshold (default:
0.75)")
parser.add_argument("--output-model", type=str, default="kidney_gate_model.keras", help="Output .keras
file path")
args = parser.parse_args()

train_path = args.train_dir or os.path.join(args.dataset_dir, "train")
val_path = args.val_dir or os.path.join(args.dataset_dir, "val")

print(f"TensorFlow Version: {tf.__version__}")
print(f"Train directory: {train_path}")
print(f"Validation directory: {val_path}")
print(f"Strict Gate Threshold: {args.threshold * 100:.0f}%")

# 1. Clean & Audit Dataset
clean_and_audit_dataset(train_path)
clean_and_audit_dataset(val_path)

# 2. tf.data Pipelines with aspect ratio preservation
img_size = (args.img_size, args.img_size)
train_raw = tf.keras.utils.image_dataset_from_directory(
    train_path,
    image_size=img_size,
    crop_to_aspect_ratio=True,
    batch_size=args.batch_size,
    label_mode="binary",
    shuffle=True,
    seed=SEED
)
val_raw = tf.keras.utils.image_dataset_from_directory(
    val_path,
    image_size=img_size,
    crop_to_aspect_ratio=True,
    batch_size=args.batch_size,
    label_mode="binary",
    shuffle=False
)

class_names = train_raw.class_names
print(f"Classes detected: {class_names}")

data_augmentation = tf.keras.Sequential([
    layers.RandomFlip("horizontal"),
    layers.RandomRotation(0.08),
    layers.RandomZoom(0.08),
    layers.RandomContrast(0.12),
    layers.RandomBrightness(0.08)
], name="gate_data_augmentation")

AUTOTUNE = tf.data.AUTOTUNE
train_ds = (
    train_raw.map(lambda x, y: (data_augmentation(x, training=True), y), num_parallel_calls=AUTOTUNE)
    .shuffle(buffer_size=1000, seed=SEED)
    .prefetch(buffer_size=AUTOTUNE)
)
val_ds = val_raw.prefetch(buffer_size=AUTOTUNE)

```

```

# 3. Build Model
gate_model, base_backbone = build_gate_model(input_shape=(args.img_size, args.img_size, 3))
gate_model.summary()

# 4. Phase 1: Feature Extraction (Head Only)
print("\n--- Phase 1: Feature Extraction ---")
gate_model.compile(
    optimizer=optimizers.Adam(learning_rate=1e-3),
    loss="binary_crossentropy",
    metrics=[
        "accuracy",
        metrics.Precision(name="precision"),
        metrics.Recall(name="recall"),
        metrics.AUC(name="auc")
    ]
)

cb_p1 = [
    callbacks.EarlyStopping(monitor="val_auc", mode="max", patience=3, restore_best_weights=True,
verbose=1),
    callbacks.ReduceLROnPlateau(monitor="val_loss", factor=0.5, patience=2, min_lr=1e-6, verbose=1)
]
gate_model.fit(train_ds, validation_data=val_ds, epochs=args.epochs_p1, callbacks=cb_p1)

# 5. Phase 2: Fine-Tuning Top 35 Layers
print("\n--- Phase 2: Fine-Tuning ---")
base_backbone.trainable = True
for layer in base_backbone.layers[:-35]:
    layer.trainable = False

gate_model.compile(
    optimizer=optimizers.Adam(learning_rate=1e-5),
    loss="binary_crossentropy",
    metrics=[
        "accuracy",
        metrics.Precision(name="precision"),
        metrics.Recall(name="recall"),
        metrics.AUC(name="auc")
    ]
)

cb_p2 = [
    callbacks.EarlyStopping(monitor="val_auc", mode="max", patience=4, restore_best_weights=True,
verbose=1),
    callbacks.ModelCheckpoint(filepath=args.output_model, monitor="val_auc", mode="max",
save_best_only=True, verbose=1),
    callbacks.ReduceLROnPlateau(monitor="val_loss", factor=0.5, patience=2, min_lr=1e-7, verbose=1)
]
gate_model.fit(train_ds, validation_data=val_ds, epochs=args.epochs_p2, callbacks=cb_p2)

# 6. Save final model
gate_model.save(args.output_model)
print(f"\n[SUCCESS] Model successfully trained and saved to: {os.path.abspath(args.output_model)}")

# 7. Verification reload
print("Verifying model reload...")
loaded = tf.keras.models.load_model(args.output_model)
dummy = np.zeros((1, args.img_size, args.img_size, 3), dtype=np.float32)
pred = loaded.predict(dummy, verbose=0)
print(f"Reload verification successful. Dummy output shape: {pred.shape}")

```

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```

"""
KidneyVision AI – Stage 2 Diagnostic Model Training Script
ResNet50V2 Transfer Learning for Nephrolithiasis (Kidney Stone) Detection with Grad-CAM

Usage:
    python train_diagnostic_model.py --dataset-dir /path/to/dataset_diagnostic --output-model
    kidney_stone_diagnostic_model.keras
"""

import os
import sys
import glob
import hashlib
import random
import argparse
import numpy as np
from PIL import Image

import tensorflow as tf
from tensorflow.keras import layers, models, callbacks, optimizers, metrics

# Reproducibility
SEED = 42
os.environ["PYTHONHASHSEED"] = str(SEED)
random.seed(SEED)
np.random.seed(SEED)
tf.random.set_seed(SEED)

def clean_and_audit_dataset(directory_path: str):
    """Scans directory to remove unreadable, zero-byte, and duplicate images."""
    if not os.path.exists(directory_path):
        print(f"[WARNING] Directory not found: {directory_path}. Skipping audit.")
        return

    valid_extensions = (".jpg", ".jpeg", ".png", ".bmp", ".tif", ".tiff", ".jfif")
    seen_hashes = set()
    removed_corrupted = 0
    removed_duplicates = 0
    total_valid = 0

    print(f"\n[DATA AUDIT] Auditing images in: {directory_path}...")

    for root, _, files in os.walk(directory_path):
        for fname in files:
            fpath = os.path.join(root, fname)

            if not fname.lower().endswith(valid_extensions):
                continue

            # Check file size > 0
            if os.path.getsize(fpath) == 0:
                try:
                    os.remove(fpath)
                    removed_corrupted += 1
                except OSError:
                    pass
                continue

            # Verify image format integrity
            try:

```

```

        with Image.open(fpath) as img:
            img.verify()

        # MD5 Hash check for duplicate removal
        with open(fpath, "rb") as f:
            file_hash = hashlib.md5(f.read()).hexdigest()

        if file_hash in seen_hashes:
            try:
                os.remove(fpath)
                removed_duplicates += 1
            except OSError:
                pass
        else:
            seen_hashes.add(file_hash)
            total_valid += 1

    except Exception:
        try:
            os.remove(fpath)
            removed_corrupted += 1
        except OSError:
            pass

print(f"Audit Complete for {directory_path}:")
print(f" - Valid unique images: {total_valid}")
print(f" - Corrupted files removed: {removed_corrupted}")
print(f" - Duplicate files removed: {removed_duplicates}")

def build_diagnostic_model(input_shape=(224, 224, 3)):
    """Builds ResNet50V2 diagnostic model optimized for Grad-CAM explainability."""
    base_model = tf.keras.applications.ResNet50V2(
        input_shape=input_shape,
        include_top=False,
        weights="imagenet"
    )
    base_model.trainable = False # Frozen for Phase 1

    inputs = tf.keras.Input(shape=input_shape, name="scan_input")
    x = tf.keras.applications.resnet_v2.preprocess_input(inputs)
    x = base_model(x, training=False)

    # Clinical Diagnostic Head with GAP
    x = layers.GlobalAveragePooling2D(name="gap")(x)
    x = layers.BatchNormalization(name="diag_bn1")(x)
    x = layers.Dropout(0.4, name="diag_dropout1")(x)
    x = layers.Dense(128, activation="relu", name="diag_dense1")(x)
    x = layers.BatchNormalization(name="diag_bn2")(x)
    x = layers.Dropout(0.2, name="diag_dropout2")(x)
    outputs = layers.Dense(1, activation="sigmoid", name="stone_probability")(x)

    model = tf.keras.Model(inputs=inputs, outputs=outputs, name="KidneyVision_ResNet50V2_Diagnostic")
    return model, base_model

def main():
    parser = argparse.ArgumentParser(description="Train KidneyVision AI Diagnostic Model")
    parser.add_argument("--dataset-dir", type=str, default="dataset_diagnostic", help="Path to diagnostic dataset directory")
    parser.add_argument("--train-dir", type=str, default=None, help="Explicit path to train directory")

```

```

parser.add_argument("--val-dir", type=str, default=None, help="Explicit path to val directory")
parser.add_argument("--img-size", type=int, default=224, help="Target image dimension (default: 224)")
parser.add_argument("--batch-size", type=int, default=32, help="Batch size (default: 32)")
parser.add_argument("--epochs-p1", type=int, default=8, help="Phase 1 epochs (frozen backbone)")
parser.add_argument("--epochs-p2", type=int, default=15, help="Phase 2 epochs (fine-tuning)")
parser.add_argument("--output-model", type=str, default="kidney_stone_diagnostic_model.keras",
help="Output .keras file path")
args = parser.parse_args()

train_path = args.train_dir or os.path.join(args.dataset_dir, "train")
val_path = args.val_dir or os.path.join(args.dataset_dir, "val")

print(f"TensorFlow Version: {tf.__version__}")
print(f"Train directory: {train_path}")
print(f"Validation directory: {val_path}")

# 1. Clean & Audit Dataset
clean_and_audit_dataset(train_path)
clean_and_audit_dataset(val_path)

# 2. tf.data Pipelines
img_size = (args.img_size, args.img_size)
train_raw = tf.keras.utils.image_dataset_from_directory(
    train_path,
    image_size=img_size,
    batch_size=args.batch_size,
    label_mode="binary",
    shuffle=True,
    seed=SEED
)
val_raw = tf.keras.utils.image_dataset_from_directory(
    val_path,
    image_size=img_size,
    batch_size=args.batch_size,
    label_mode="binary",
    shuffle=False
)

class_names = train_raw.class_names
print(f"Classes detected: {class_names}")

data_augmentation = tf.keras.Sequential([
    layers.RandomFlip("horizontal"),
    layers.RandomRotation(0.06),
    layers.RandomZoom(0.06),
    layers.RandomContrast(0.12),
    layers.RandomBrightness(0.08)
], name="medical_augmentation")

AUTOTUNE = tf.data.AUTOTUNE
train_ds = (
    train_raw.map(lambda x, y: (data_augmentation(x, training=True), y), num_parallel_calls=AUTOTUNE)
    .shuffle(buffer_size=1000, seed=SEED)
    .prefetch(buffer_size=AUTOTUNE)
)
val_ds = val_raw.prefetch(buffer_size=AUTOTUNE)

# 3. Build Model
diag_model, base_backbone = build_diagnostic_model(input_shape=(args.img_size, args.img_size, 3))
diag_model.summary()

```

```

# 4. Phase 1: Feature Extraction (Diagnostic Head)
print("\n--- Phase 1: Feature Extraction ---")
diag_model.compile(
    optimizer=optimizers.Adam(learning_rate=1e-3),
    loss="binary_crossentropy",
    metrics=[
        "accuracy",
        metrics.Precision(name="precision"),
        metrics.Recall(name="recall"),
        metrics.AUC(name="auc")
    ]
)

cb_p1 = [
    callbacks.EarlyStopping(monitor="val_auc", mode="max", patience=3, restore_best_weights=True,
verbose=1),
    callbacks.ReduceLROnPlateau(monitor="val_loss", factor=0.5, patience=2, min_lr=1e-6, verbose=1)
]
diag_model.fit(train_ds, validation_data=val_ds, epochs=args.epochs_p1, callbacks=cb_p1)

# 5. Phase 2: Fine-Tuning Top Residual Layers
print("\n--- Phase 2: Fine-Tuning ---")
base_backbone.trainable = True
for layer in base_backbone.layers[:-40]:
    layer.trainable = False

diag_model.compile(
    optimizer=optimizers.Adam(learning_rate=1e-5),
    loss="binary_crossentropy",
    metrics=[
        "accuracy",
        metrics.Precision(name="precision"),
        metrics.Recall(name="recall"),
        metrics.AUC(name="auc")
    ]
)

cb_p2 = [
    callbacks.EarlyStopping(monitor="val_auc", mode="max", patience=4, restore_best_weights=True,
verbose=1),
    callbacks.ModelCheckpoint(filepath=args.output_model, monitor="val_auc", mode="max",
save_best_only=True, verbose=1),
    callbacks.ReduceLROnPlateau(monitor="val_loss", factor=0.5, patience=2, min_lr=1e-7, verbose=1)
]
diag_model.fit(train_ds, validation_data=val_ds, epochs=args.epochs_p2, callbacks=cb_p2)

# 6. Save final model
diag_model.save(args.output_model)
print(f"\n[SUCCESS] Model successfully trained and saved to: {os.path.abspath(args.output_model)}")

# 7. Verification reload
print("Verifying model reload...")
loaded = tf.keras.models.load_model(args.output_model)
dummy = np.zeros((1, args.img_size, args.img_size, 3), dtype=np.float32)
pred = loaded.predict(dummy, verbose=0)
print(f"Reload verification successful. Dummy output shape: {pred.shape}")

if __name__ == "__main__":
    main()

```

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