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```

def is_port_in_use(port: int, host: str = "127.0.0.1") -> bool:
    with socket.socket(socket.AF_INET, socket.SOCK_STREAM) as s:
        s.settimeout(0.5)
        return s.connect_ex((host, port)) == 0

def kill_pids_on_ports(ports):
    """Finds and kills any processes holding specified ports on Windows."""
    if sys.platform != "win32":
        return
    for port in ports:
        try:
            cmd = f'netstat -ano | findstr /C:"{port}" '
            output = subprocess.check_output(cmd, shell=True, text=True, stderr=subprocess.DEVNULL)
            for line in output.splitlines():
                if "LISTENING" in line:
                    parts = line.strip().split()
                    pid = parts[-1]
                    if pid and pid != "0":
                        subprocess.run(
                            ["taskkill", "/F", "/PID", pid],
                            stdout=subprocess.DEVNULL,
                            stderr=subprocess.DEVNULL,
                        )
                        log("RUNNER", f"Freed port {port} (killed PID {pid})", YELLOW)
        except Exception:
            pass

def find_python_for_ai() -> str:
    # Prefer virtual environment in ai-flask-service
    venv_win = FLASK_DIR / "venv" / "Scripts" / "python.exe"
    venv_unix = FLASK_DIR / "venv" / "bin" / "python"
    if venv_win.exists():
        return str(venv_win)
    if venv_unix.exists():
        return str(venv_unix)
    # Check root venv
    root_win = ROOT_DIR / "venv" / "Scripts" / "python.exe"
    if root_win.exists():
        return str(root_win)
    return sys.executable

def find_php() -> str:
    # 1. Check current PATH
    php = shutil.which("php")
    if php:
        return php

    # 2. Check WinGet package location on Windows
    local_app_data = os.environ.get("LOCALAPPDATA", "")
    if local_app_data:
        winget_pattern = Path(local_app_data) / "Microsoft" / "WinGet" / "Packages"
        if winget_pattern.exists():
            matches = list(winget_pattern.glob("PHP.PHP*/**/php.exe"))
            if matches and matches[0].exists():
                php_found = str(matches[0])
                php_dir = str(matches[0].parent)
                if php_dir not in os.environ.get("PATH", ""):
                    os.environ["PATH"] = f"{php_dir};{os.environ.get('PATH', '')}"

```

```

        return php_found

# 3. Common Laragon or XAMPP fallback locations on Windows
common_paths = [
    Path(r"C:\laragon\bin\php\php-8.3.30-Win32-vs16-x64\php.exe"),
    Path(r"C:\laragon\bin\php\php-8.2.0-Win32-vs16-x64\php.exe"),
    Path(r"C:\xampp\php\php.exe"),
]
for p in common_paths:
    if p.exists():
        php_dir = str(p.parent)
        if php_dir not in os.environ.get("PATH", ""):
            os.environ["PATH"] = f"{php_dir};{os.environ.get('PATH', '')}"
        return str(p)

return "php"

def find_npm() -> str:
    npm = shutil.which("npm.cmd") if sys.platform == "win32" else shutil.which("npm")
    return npm if npm else "npm"

def setup_environments():
    log("RUNNER", "Checking configuration & environment files...", BOLD)

    # 1. Root .env
    root_env = ROOT_DIR / ".env"
    root_example = ROOT_DIR / ".env.example"
    if not root_env.exists() and root_example.exists():
        shutil.copy(root_example, root_env)
        log("RUNNER", "Created root .env from .env.example", GREEN)

    # 2. Flask .env
    flask_env = FLASK_DIR / ".env"
    flask_example = FLASK_DIR / ".env.example"
    if not flask_env.exists() and flask_example.exists():
        shutil.copy(flask_example, flask_env)
        log("RUNNER", "Created ai-flask-service .env", GREEN)

    # 3. Laravel .env & SQLite DB
    laravel_env = LARAVEL_DIR / ".env"
    laravel_example = LARAVEL_DIR / ".env.example"
    if not laravel_env.exists() and laravel_example.exists():
        shutil.copy(laravel_example, laravel_env)
        log("RUNNER", "Created backend-laravel .env", GREEN)

    db_file = LARAVEL_DIR / "database" / "database.sqlite"
    if not db_file.exists():
        db_file.parent.mkdir(parents=True, exist_ok=True)
        db_file.touch()
        log("RUNNER", f"Initialized SQLite database at {db_file}", GREEN)

    # 4. Frontend .env
    frontend_env = FRONTEND_DIR / ".env"
    frontend_example = FRONTEND_DIR / ".env.example"
    if not frontend_env.exists() and frontend_example.exists():
        shutil.copy(frontend_example, frontend_env)
        log("RUNNER", "Created frontend-react .env", GREEN)

```

```

def stream_output(pipe, prefix: str, color: str):
    try:
        for line in iter(pipe.readline, ""):
            if not line:
                break
            line_str = line.strip()
            if line_str:
                log(prefix, line_str, color)
    except Exception:
        pass
    finally:
        pipe.close()

def kill_proc(proc):
    try:
        if sys.platform == "win32":
            subprocess.run(
                ["taskkill", "/F", "/T", "/PID", str(proc.pid)],
                stdout=subprocess.DEVNULL,
                stderr=subprocess.DEVNULL,
            )
        else:
            proc.terminate()
            try:
                proc.wait(timeout=3)
            except subprocess.TimeoutExpired:
                proc.kill()
    except Exception:
        pass

def shutdown_all(*args):
    global SHUTTING_DOWN
    if SHUTTING_DOWN:
        return
    SHUTTING_DOWN = True
    print("\n")
    log("RUNNER", "Stopping all KidneyVision AI microservices...", YELLOW)
    for p in PROCESSES:
        kill_proc(p)
    log("RUNNER", "All services stopped cleanly. Goodbye!", GREEN)
    sys.exit(0)

def wait_for_service(url: str, timeout: int = 45) -> bool:
    import urllib.request
    start_time = time.time()
    while time.time() - start_time < timeout:
        if SHUTTING_DOWN:
            return False
        try:
            req = urllib.request.Request(url, headers={"User-Agent": "KidneyVisionRunner/1.0"})
            with urllib.request.urlopen(req, timeout=2) as response:
                if response.status in (200, 301, 302, 304):
                    return True
        except Exception:
            time.sleep(1)
    return False

```

```

def monitor_readiness(open_browser: bool):
    time.sleep(2)
    log("RUNNER", "Waiting for services to become available...", CYAN)

    ai_ready = wait_for_service("http://127.0.0.1:5000/health", timeout=35)
    backend_ready = wait_for_service("http://127.0.0.1:8000/api/health", timeout=35)
    frontend_ready = wait_for_service("http://localhost:3000", timeout=35)

    if SHUTTING_DOWN:
        return

    banner = (
        f"\n{BOLD}{GREEN}===== {RESET}\n"
        f"{BOLD}          KIDNEYVISION AI - ALL SERVICES ARE RUNNING! {RESET}\n"
        f"{BOLD}{GREEN}===== {RESET}\n"
        f" {CYAN}* React Clinical Dashboard: {RESET} {BOLD}http://localhost:3000 {RESET}\n"
        f" {MAGENTA}* Laravel Backend API: {RESET} {BOLD}http://127.0.0.1:8000/api {RESET}\n"
        f" {YELLOW}* Flask AI Inference Engine: {RESET} {BOLD}http://127.0.0.1:5000/health {RESET}\n"
        f"{BOLD}{GREEN}===== {RESET}\n"
        f" Press {BOLD}Ctrl + C {RESET} in this terminal anytime to cleanly shut down all services.\n"
    )
    print(banner, flush=True)

    if open_browser and frontend_ready:
        try:
            webbrowser.open("http://localhost:3000")
        except Exception:
            pass

def main():
    parser = argparse.ArgumentParser(description="Run all KidneyVision AI microservices.")
    parser.add_argument(
        "--no-browser",
        action="store_true",
        help="Do not open the web browser automatically",
    )
    parser.add_argument(
        "--kill-existing",
        action="store_true",
        help="Automatically terminate any existing processes using ports 3000, 5000, or 8000",
    )
    args = parser.parse_args()

    # Register exit handlers
    signal.signal(signal.SIGINT, shutdown_all)
    signal.signal(signal.SIGTERM, shutdown_all)

    print(f"\n{BOLD}{CYAN}=== KidneyVision AI Platform Orchestrator === {RESET}\n")

    setup_environments()

    ports_map = {
        5000: "AI Flask Microservice",
        8000: "Laravel Backend API",
        3000: "React Frontend",
    }

    if args.kill_existing:
        kill_pids_on_ports(list(ports_map.keys()))

```

```

# Port availability check
for port, name in ports_map.items():
    if is_port_in_use(port):
        log(
            "WARN",
            f"Port {port} ({name}) is in use! If a conflict occurs, pass --kill-existing to free it.",
            YELLOW,
        )

ai_python = find_python_for_ai()
php_bin = find_php()
npm_bin = find_npm()

log("RUNNER", f"Python (AI): {ai_python}", CYAN)
log("RUNNER", f"PHP (Laravel): {php_bin}", CYAN)
log("RUNNER", f"NPM (Frontend): {npm_bin}", CYAN)

# 1. Start AI Flask Service
log("RUNNER", "Starting AI Flask Microservice on http://127.0.0.1:5000 ...", YELLOW)
proc_flask = subprocess.Popen(
    [ai_python, "app.py"],
    cwd=str(FLASK_DIR),
    stdout=subprocess.PIPE,
    stderr=subprocess.STDOUT,
    text=True,
    bufsize=1,
)
PROCESSES.append(proc_flask)
threading.Thread(
    target=stream_output,
    args=(proc_flask.stdout, "AI-FLASK", YELLOW),
    daemon=True,
).start()

# 2. Start Laravel Backend API
log("RUNNER", "Starting Laravel Backend API on http://127.0.0.1:8000 ...", MAGENTA)
proc_laravel = subprocess.Popen(
    [php_bin, "artisan", "serve", "--host=127.0.0.1", "--port=8000"],
    cwd=str(LARAVEL_DIR),
    stdout=subprocess.PIPE,
    stderr=subprocess.STDOUT,
    text=True,
    bufsize=1,
)
PROCESSES.append(proc_laravel)
threading.Thread(
    target=stream_output,
    args=(proc_laravel.stdout, "BACKEND", MAGENTA),
    daemon=True,
).start()

# 3. Start React Frontend Client
log("RUNNER", "Starting React Frontend Client on http://localhost:3000 ...", CYAN)
proc_frontend = subprocess.Popen(
    [npm_bin, "run", "dev"],
    cwd=str(FRONTEND_DIR),
    stdout=subprocess.PIPE,
    stderr=subprocess.STDOUT,
    text=True,
    bufsize=1,
)

```

```

PROCESSES.append(proc_frontend)
threading.Thread(
    target=stream_output,
    args=(proc_frontend.stdout, "FRONTEND", CYAN),
    daemon=True,
).start()

# 4. Start health check and banner monitor
monitor_thread = threading.Thread(
    target=monitor_readiness,
    args=(not args.no_browser,),
    daemon=True,
)
monitor_thread.start()

# Keep main thread alive and watch child processes
process_names = {
    proc_flask: "AI Flask Microservice",
    proc_laravel: "Laravel Backend API",
    proc_frontend: "React Frontend Client",
}
reported_exits = set()
try:
    while True:
        for p, name in process_names.items():
            code = p.poll()
            if code is not None and not SHUTTING_DOWN and p not in reported_exits:
                reported_exits.add(p)
                log("WARN", f"[{name}] exited unexpectedly with code {code}.", RED)
                time.sleep(1)
except KeyboardInterrupt:
    shutdown_all()

if __name__ == "__main__":
    main()

```

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```

version: '3.8'

services:
  # Flask AI Deep Learning Microservice
  ai-service:
    build:
      context: .
      dockerfile: deployment/docker/Dockerfile.ai-service
    container_name: kidneyvision-ai-service
    restart: unless-stopped
    ports:
      - "5000:5000"
    environment:
      - FLASK_ENV=production
      - MODEL_PATH=best_model.keras
      - PORT=5000
    volumes:
      - ./models/production:/app/models
    healthcheck:
      test: ["CMD", "curl", "-f", "http://localhost:5000/health"]
      interval: 30s
      timeout: 10s
      retries: 3

  # Laravel 12 Backend API
  backend:
    build:
      context: .
      dockerfile: deployment/docker/Dockerfile.backend
    container_name: kidneyvision-backend
    restart: unless-stopped
    ports:
      - "8000:8000"
    environment:
      - APP_ENV=production
      - APP_DEBUG=false
      - APP_KEY=${APP_KEY:-base64:cK15Y1E4Y00F0GQd3H0Y5B40pZ1vB5h0+H1f1m9X+8M=}
      - APP_URL=http://localhost:8000
      - FRONTEND_URL=http://localhost:3000
      - DB_CONNECTION=sqlite
      - DB_DATABASE=/var/www/html/database/database.sqlite
      - AI_SERVICE_URL=http://ai-service:5000
    depends_on:
      ai-service:
        condition: service_healthy
    command: php artisan serve --host=0.0.0.0 --port=8000
    volumes:
      - ./backend-laravel/storage:/var/www/html/storage
      - ./backend-laravel/database:/var/www/html/database

  # React 19 Frontend Web Client
  frontend:
    build:
      context: .
      dockerfile: deployment/docker/Dockerfile.frontend
    container_name: kidneyvision-frontend
    restart: unless-stopped

```

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